

Anti-SLITRK6 antibody [SAA0827] (STJA0027697)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.14 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	30 mM Histidine, pH5.8, 10% sucrose, 0.02% Tween80
Isotype	IgG2k
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	84189
Gene Symbol	SLITRK6
UniProt ID	SLIK6_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--