

Anti-CD26/DPP4 antibody [SAA1395] (STJA0027727)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 1 mg/ml |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant. |
| Formulation | 0.01M PBS, pH 7.4, 0.09% Sodium azide. |
| Isotype | IgG2ak |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

TARGET INFORMATION

| | |
|-------------|--------------------------|
| Gene ID | 25253 |
| Gene Symbol | Dpp4 |
| UniProt ID | DPP4_RAT |

ADDITIONAL INFORMATION

| | |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|