

**Anti-CSPG4/HMW-MAA/MCSP/NG2 antibody [763.74#] (STJA0027780)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1.22 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">1464</a>
Gene Symbol	<a href="#">CSPG4</a>
UniProt ID	<a href="#">CSPG4_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--