

## Anti-CD121a/IL1R1 antibody [15C4] (STJA0027805)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/IHC/SPR/WB
Host / Source	Mouse
Reactivity	Human/Monkey

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Isotype	IgG2ak
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">3554</a>
Gene Symbol	<a href="#">IL1R1</a>
UniProt ID	<a href="#">IL1R1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--