

Anti-Complement C5b-9 Complex antibody [aE11#] (STJA0027851)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IHC/WB
Host / Source	Mouse
Reactivity	Human/Baboon/Horse/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.93 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2ak
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--