

## Anti-PRLR antibody [SAA2060] (STJA0027854)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FC           |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                          |
| Concentration       | 1 mg/ml                                                                                             |
| Conjugation         | Unconjugated                                                                                        |
| Purification        | Protein A/G purified from cell culture supernatant.                                                 |
| Formulation         | 0.01M PBS, pH 7.4, 0.09% Sodium azide.                                                              |
| Isotype             | IgG2ak                                                                                              |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">5618</a>       |
| Gene Symbol | <a href="#">PRLR</a>       |
| UniProt ID  | <a href="#">PRLR_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|