

Anti-NECTIN4/PVRL4 antibody [Nec4.05] (STJA0027865)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC |
| Host / Source | Human |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 1 mg/ml |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant. |
| Formulation | 0.01M PBS, pH 7.4, 0.09% Sodium azide. |
| Isotype | IgG1k |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 81607 |
| Gene Symbol | NECTIN4 |
| UniProt ID | NECT4_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|