

Anti-CDw210b/IL10RB antibody [SAA2106] (STJA0027875)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Isotype	VHH-mFc (IgG2a)
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	3588
Gene Symbol	IL10RB
UniProt ID	I10R2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--