

Anti-Mycobacterium tuberculosis (strain ATCC 25618/H37Rv) antibody [SAA2170] (STJA0027890)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/SPR
Host / Source	Human
Reactivity	M.tuberculosis

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--