

**Anti-CD156c/ADAM10 antibody [8C7] {FITC} (STJA0027917)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Bovine

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

**TARGET INFORMATION**

Gene ID	<a href="#">11487</a>
Gene Symbol	<a href="#">Adam10</a>
UniProt ID	<a href="#">ADA10_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--