

**Anti-CD262/TNFRSF10B/TRAIL-R2 antibody [SAA1413] {FITC} (STJA0027920)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

**TARGET INFORMATION**

Gene ID	<a href="#">8795</a>
Gene Symbol	<a href="#">TNFRSF10B</a>
UniProt ID	<a href="#">TR10B_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--