

**Anti-CD8a/Lyt2 antibody [19178#] {FITC} (STJA0027998)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Rat
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

**TARGET INFORMATION**

Gene ID	<a href="#">12525</a>
Gene Symbol	<a href="#">Cd8a</a>
UniProt ID	<a href="#">CD8A_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--