

**Anti-CD18/ITGB2 antibody [SAA0815] {FITC} (STJA0028038)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Human
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

**TARGET INFORMATION**

Gene ID	<a href="#">3689</a>
Gene Symbol	<a href="#">ITGB2</a>
UniProt ID	<a href="#">ITB2_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------