

## Anti-CD2 antibody [SAA0001] {FITC} (STJA0028062)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

### TARGET INFORMATION

Gene ID	<a href="#">914</a>
Gene Symbol	<a href="#">CD2</a>
UniProt ID	<a href="#">CD2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--