

## Anti-LPL/LIPD antibody [5D2] {FITC} (STJA0028073)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF-P/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Bovine

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

### TARGET INFORMATION

Gene ID	<a href="#">4023</a>
Gene Symbol	<a href="#">LPL</a>
UniProt ID	<a href="#">LIPL_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------