

Anti-OPN/SPP1 antibody [SAA0512] {FITC} (STJA0028151)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/FC/WB
	Host / Source	Mouse
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	05 mg/mL
	Conjugation	FITC
	Purification	Purified
	Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
	Isotype	IgG2ak
	Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze
TARGET INFORMATION		
	Gene ID	6696
	Gene Symbol	SPP1
	UniProt ID	OSTP_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS