

Anti-CD11b/ITGAM antibody [DC13] {FITC} (STJA0028166)**GENERAL INFORMATION**

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC/IF |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 05 mg/mL |
| Conjugation | FITC |
| Purification | Purified |
| Formulation | PBS pH 7.4, 50% glycerol, 0.05% proclin300 |
| Isotype | VHH-rFc (IgG2b) |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

TARGET INFORMATION

| | |
|-------------|--|
| Gene Symbol | Itgam Uniprot ID= ITAM_MOUSE " |
|-------------|--|

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|