

**Anti-CD324/CDH1 antibody [SAA0269] {FITC} (STJA0028188)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | human              |

**PRODUCT PROPERTIES**

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Clonality           | Monoclonal                                                  |
| Concentration       | 05 mg/mL                                                    |
| Conjugation         | FITC                                                        |
| Purification        | Purified                                                    |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                  |
| Isotype             | IgG2ak                                                      |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">999</a>         |
| Gene Symbol | <a href="#">CDH1</a>        |
| UniProt ID  | <a href="#">CADH1_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|