

**Anti-Recombinant-CD49d/ITGA4 antibody [SAA0019] {FITC} (STJA0028194)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgGk
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

**TARGET INFORMATION**

Gene ID	<a href="#">3676</a>
Gene Symbol	<a href="#">ITGA4</a>
UniProt ID	<a href="#">ITA4_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--