

## Anti-FOLR1 antibody [SAA0118] {FITC} (STJA0028220)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Human              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 05 mg/mL  |
| Conjugation         | FITC  |
| Purification        | Purified  |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                  |
| Isotype             | IgG1k   |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2348</a>        |
| Gene Symbol | <a href="#">FOLR1</a>       |
| UniProt ID  | <a href="#">FOLR1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|