

## Anti-CD1d antibody [1B1] {FITC} (STJA0028234)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2b
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

### TARGET INFORMATION

Gene ID	<a href="#">12479</a>
Gene Symbol	<a href="#">Cd1d1</a>
UniProt ID	<a href="#">CD1D1_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------