

**Anti-CD40/TNFRSF5 antibody [FGK45] {FITC} (STJA0028364)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

**PRODUCT PROPERTIES**

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Clonality           | Monoclonal                                                  |
| Concentration       | 05 mg/mL                                                    |
| Conjugation         | FITC                                                        |
| Purification        | Purified                                                    |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                  |
| Isotype             | IgG2ak                                                      |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">21939</a>       |
| Gene Symbol | <a href="#">Cd40</a>        |
| UniProt ID  | <a href="#">TNFR5_MOUSE</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|