

Anti-CD26/DPP4 antibody [K9/9] {FITC} (STJA0028376)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/WB
Host / Source	Rat
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

TARGET INFORMATION

Gene ID	13482
Gene Symbol	Dpp4
UniProt ID	DPP4_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--