

Anti-Calreticulin/CALR antibody [SAA0511] {FITC} (STJA0028378)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

TARGET INFORMATION

Gene ID	811
Gene Symbol	CALR
UniProt ID	CALR_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--