

Anti-IL13 antibody [SAA0383] {FITC} (STJA0028427)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG4L
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

TARGET INFORMATION

Gene ID	3596
Gene Symbol	IL13
UniProt ID	IL13_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--