

Anti-CD309/KDR/VEGFR-2 antibody [SSS-23] {FITC} (STJA0028432)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	05 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

TARGET INFORMATION

Gene ID	3791
Gene Symbol	KDR
UniProt ID	VGFR2 HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--