

**Anti-GPC3 antibody [SAA0134] {FITC} (STJA0028496)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FC/WB        |
| Host / Source | Human              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Clonality           | Monoclonal                                                  |
| Concentration       | 073 mg/mL                                                   |
| Conjugation         | FITC                                                        |
| Purification        | Purified                                                    |
| Formulation         | Liquid in PBS containing 0.2% BSA and 0.05% Proclin300.     |
| Isotype             | IgG1k                                                       |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">2719</a>       |
| Gene Symbol | <a href="#">GPC3</a>       |
| UniProt ID  | <a href="#">GPC3_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|