

Anti-F4/80 antibody [SAA2157] {FITC} (STJA0028598)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	13733
Gene Symbol	Adgre1
UniProt ID	AGRE1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--