

Anti-CD360/IL21R antibody [18A5] {FITC} (STJA0028710)

GENERAL INFORMATION

| | |
|---------------|------------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Human/Mouse/Rat/Monkey |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | FITC |
| Purification | Purified |
| Formulation | PBS pH 7.4, 50% glycerol, 0.05% proclin300 |
| Isotype | IgG1L |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 50615 |
| Gene Symbol | IL21R |
| UniProt ID | IL21R_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|