

## Anti-TREM2 antibody [SAA0548] {FITC} (STJA0028733)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Mouse
Reactivity	Human/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">54209</a>
Gene Symbol	<a href="#">TREM2</a>
UniProt ID	<a href="#">TREM2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--