

**Anti-CCL24/Eotaxin-2/MPIF-2 antibody [SAA0461] {APC} (STJA0028871)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FC           |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.5 mg/mL  |
| Conjugation         | APC  |
| Purification        | Purified   |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                     |
| Isotype             | IgG2ak   |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">6369</a>        |
| Gene Symbol | <a href="#">CCL24</a>       |
| UniProt ID  | <a href="#">CCL24_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|