

Anti-CD156c/ADAM10 antibody [8C7] {APC} (STJA0028883)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Bovine

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	APC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	11487
Gene Symbol	Adam10
UniProt ID	ADA10_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--