

## Anti-LPL/LIPD antibody [5D2] {APC} (STJA0029039)

### GENERAL INFORMATION

|               |                        |
|---------------|------------------------|
| Product Type  | Primary antibodies     |
| Applications  | ELISA/FC/IF-P/WB       |
| Host / Source | Mouse                  |
| Reactivity    | Human/Mouse/Rat/Bovine |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.5 mg/mL  |
| Conjugation         | APC  |
| Purification        | Purified   |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                     |
| Isotype             | IgG2a  |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4023</a>       |
| Gene Symbol | <a href="#">LPL</a>        |
| UniProt ID  | <a href="#">LIPL_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|