

Anti-CD121a/IL1R1 antibody [15C4] {APC} (STJA0029175)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/WB
Host / Source	Mouse
Reactivity	Human/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	APC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	3554
Gene Symbol	IL1R1
UniProt ID	IL1R1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--