

## Anti-LGALS3 antibody [P3] {APC} (STJA0029244)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/WB
Host / Source	Rat
Reactivity	Human/Mouse/Rabbit

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	APC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

### TARGET INFORMATION

Gene Symbol	<a href="#">Lgals3</a> Uniprot ID= <a href="#">LEG3_MOUSE</a> "
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--