

**Anti-CD21/CD35 antibody [7G6] {APC} (STJA0029260)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Mouse
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	APC
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene Symbol	<a href="#">Cr2</a> Uniprot ID= <a href="#">CR2_MOUSE</a> "
-------------	---

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--