

**Anti-CD183/CXCR3 antibody [7H5] {APC} (STJA0029455)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                |
|---------------------|----------------------------------------------------------------|
| Clonality           | Monoclonal                                                     |
| Concentration       | 0.5 mg/mL                                                      |
| Conjugation         | APC                                                            |
| Purification        | Purified                                                       |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                     |
| Isotype             | IgG2ak                                                         |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2833</a>        |
| Gene Symbol | <a href="#">CXCR3</a>       |
| UniProt ID  | <a href="#">CXCR3_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|