

## Anti-CD200R1/OX2R antibody [OX-108] {APC} (STJA0029631)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.5 mg/mL  |
| Conjugation         | APC  |
| Purification        | Purified   |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                     |
| Isotype             | IgG1k  |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">131450</a>      |
| Gene Symbol | <a href="#">CD200R1</a>     |
| UniProt ID  | <a href="#">MO2R1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|