

**Anti-SLC34A2/NaPi2b antibody [SAA0766] {PerCP} (STJA0029892)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/IF/WB
Host / Source	Human
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">10568</a>
Gene Symbol	<a href="#">SLC34A2</a>
UniProt ID	<a href="#">NPT2B_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--