

Anti-CD230 antibody [POM1] {PerCP} (STJA0029954)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/FC
	Host / Source	Mouse
	Reactivity	Human/Mouse/Bovine/Hamster
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	0.5 mg/mL
	Conjugation	PerCP
	Purification	Purified
	Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
	Isotype	IgG2ak
	Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.
TARGET INFORMATION		
	Gene ID	19122
	Gene Symbol	Prnp
	UniProt ID	PRIO_MOUSE
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.