

## Anti-CD28 antibody [SAA0010] {PerCP} (STJA0030089)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG4
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">940</a>
Gene Symbol	<a href="#">CD28</a>
UniProt ID	<a href="#">CD28_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------