

Anti-CD204/MSR1 antibody [H1H21234] {PerCP} (STJA0030251)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	4481
Gene Symbol	MSR1
UniProt ID	MSRE_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--