

**Anti-CD32b/FCGR2B antibody [SAA0547] {PerCP} (STJA0030339)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Human
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">2213</a>
Gene Symbol	<a href="#">FCGR2B</a>
UniProt ID	<a href="#">FCG2B_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--