

Anti-CDw210b/IL10RB antibody [SAA2106] {PerCP} (STJA0030496)**GENERAL INFORMATION**

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | PerCP |
| Purification | Purified |
| Formulation | PBS pH 7.4, 50% glycerol, 0.05% proclin300 |
| Isotype | VHH-mFc (IgG2a) |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 3588 |
| Gene Symbol | IL10RB |
| UniProt ID | I10R2_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|