

Anti-CD279/PDCD1/PD1 antibody [SAA1476] {PerCP} (STJA0030541)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Dog
Reactivity	Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--