

**Anti-GOLM1/GP73 antibody [SAA0106] {PerCP} (STJA0030591)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">51280</a>
Gene Symbol	<a href="#">GOLM1</a>
UniProt ID	<a href="#">GOLM1_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--