

Anti-CD2/LFA-2/OX-34 antibody [OX-34] {PE} (STJA0030961)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IF/WB
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PE
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	497761
Gene Symbol	Cd2
UniProt ID	CD2_RAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--