

## Anti-Recombinant-CD62L/SELL antibody [SAA1992] {PE} (STJA0031100)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.5 mg/mL  |
| Conjugation         | PE   |
| Purification        | Purified   |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                     |
| Isotype             | IgG  |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">6402</a>        |
| Gene Symbol | <a href="#">SELL</a>        |
| UniProt ID  | <a href="#">LYAM1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|