

## Anti-CD36 antibody [JL2] {PE} (STJA0031165)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FC/WB        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                |
|---------------------|----------------------------------------------------------------|
| Clonality           | Monoclonal                                                     |
| Concentration       | 0.5 mg/mL                                                      |
| Conjugation         | PE                                                             |
| Purification        | Purified                                                       |
| Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                     |
| Isotype             | IgG1k                                                          |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">948</a>        |
| Gene Symbol | <a href="#">CD36</a>       |
| UniProt ID  | <a href="#">CD36_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|