

Anti-CD49f/ITGA6 antibody [P2D10] {PE} (STJA0031237)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | PE |
| Purification | Purified |
| Formulation | PBS pH 7.4, 50% glycerol, 0.05% proclin300 |
| Isotype | IgG1k |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 3655 |
| Gene Symbol | ITGA6 |
| UniProt ID | ITA6_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|